

Product Datasheet

Mouse bACE1 ELISA Kit (orb1201229)

Description	Mouse Beta-Site APP Cleaving Enzyme 1 (bACE1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid or other biological fluids.
Uniprot ID	P56818
Sensitivity	5.9pg/mL
Entrez	23821
Expiration Date	Please enquire.