

## Product Datasheet

### Human TLR2 ELISA Kit (orb1201312)

<b>Description</b>	Human Toll Like Receptor 2 (TLR2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, saliva, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O60603</a>
<b>Sensitivity</b>	0.113ng/mL
<b>Entrez</b>	<a href="#">7097</a>
<b>Expiration Date</b>	Please enquire.