

Product Datasheet

Human TFR2 ELISA Kit (orb1201319)

Description	Human Transferrin Receptor 2 (TFR2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UP52
Sensitivity	0.07ng/mL
Entrez	7036
Expiration Date	Please enquire.