

Product Datasheet

Human PDCD1LG1 ELISA Kit (orb1201409)

Description	Human Programmed Cell Death Protein 1 Ligand 1 (PDCD1LG1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZQ7
Sensitivity	0.066ng/mL
Entrez	29126
Expiration Date	Please enquire.