

## Product Datasheet

### Human MTL ELISA Kit (orb1201456)

<b>Description</b>	Human Motilin (MTL) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P12872</a>
<b>Sensitivity</b>	4.28pg/mL
<b>Entrez</b>	<a href="#">4295</a>
<b>Expiration Date</b>	Please enquire.