

## Product Datasheet

### Human MCP2 ELISA Kit (orb1201481)

<b>Description</b>	Human Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P80075</b>
<b>Sensitivity</b>	8.6pg/mL
<b>Entrez</b>	<b>6355</b>
<b>Expiration Date</b>	Please enquire.