

Product Datasheet

Human KLK6 ELISA Kit (orb1201511)

Description	Human Kallikrein 6 (KLK6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid, breast milk or other biological fluids.
Uniprot ID	Q92876
Sensitivity	0.126ng/mL
Entrez	5653
Expiration Date	Please enquire.