

Product Datasheet

Human CCR2 ELISA Kit (orb1201677)

Description	Human Chemokine C-C-Motif Receptor 2 (CCR2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P41597
Sensitivity	0.061ng/mL
Entrez	729230
Expiration Date	Please enquire.