

Product Datasheet

Human CCL3L1 ELISA Kit (orb1201680)

Description	Human Chemokine C-C-Motif Ligand 3 Like Protein 1 (CCL3L1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P16619
Sensitivity	32pg/mL
Entrez	6349
Expiration Date	Please enquire.