

Product Datasheet

Rat IL33 ELISA Kit (orb1201767)

Description	Rat Interleukin 33 (IL33) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q66H70
Sensitivity	27pg/mL
Entrez	361749
Expiration Date	Please enquire.