

Product Datasheet

Canine TNFb ELISA Kit (orb1201785)

Description	Canine Tumor Necrosis Factor Beta (TNFb) ELISA Kit
Reactivity	Canine
Range	15.6-1000ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5WR07
Sensitivity	3.6ng/mL
Entrez	607183
Expiration Date	Please enquire.