

## Product Datasheet

### Canine TGFb2 ELISA Kit (orb1201787)

<b>Description</b>	Canine Transforming Growth Factor Beta 2 (TGFb2) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, cell Lysates, cell culture supernatants and other biological fluids
<b>Uniprot ID</b>	<a href="#">F1PKH2</a>
<b>Sensitivity</b>	5.4pg/mL
<b>Entrez</b>	<a href="#">488596</a>
<b>Expiration Date</b>	Please enquire.