

Product Datasheet

Canine TGFb1 ELISA Kit (orb1201789)

Description	Canine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids
Uniprot ID	P54831
Sensitivity	5.6pg/mL
Entrez	403998
Expiration Date	Please enquire.