

Product Datasheet

Guinea IL4 ELISA Kit (orb1201801)

Description	Guinea pig Interleukin 4 (IL4) ELISA Kit
Reactivity	Guinea pig
Range	1.56-100pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	G0WM66
Sensitivity	0.57pg/mL
Entrez	100720403
Expiration Date	Please enquire.