

Product Datasheet

Canine IL18 ELISA Kit (orb1201814)

Description	Canine Interleukin 18 (IL18) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XSR0
Sensitivity	5.9pg/mL
Entrez	403796
Expiration Date	Please enquire.