

Product Datasheet

Guinea IL10 ELISA Kit (orb1201820)

Description	Guinea pig Interleukin 10 (IL10) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	H0VDI2
Sensitivity	5.2pg/mL
Entrez	100715618
Expiration Date	Please enquire.