

Product Datasheet

Canine IFN γ ELISA Kit (orb1201825)

Description	Canine Interferon Gamma (IFN γ) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42161
Sensitivity	5.9pg/mL
Entrez	403801
Expiration Date	Please enquire.