

## Product Datasheet

### Canine IFN $\beta$ ELISA Kit (orb1201826)

<b>Description</b>	Canine Interferon Beta (IFN $\beta$ ) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">B6E116</a>
<b>Sensitivity</b>	5.6pg/mL
<b>Entrez</b>	<a href="#">481558</a>
<b>Expiration Date</b>	Please enquire.