

Product Datasheet

Canine GMCSF ELISA Kit (orb1201829)

Description	Canine Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48749
Sensitivity	7.2pg/mL
Entrez	403923
Expiration Date	Please enquire.