

Product Datasheet

Canine BDNF ELISA Kit (orb1201842)

Description	Canine Brain Derived Neurotrophic Factor (BDNF) ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q7YRB4
Sensitivity	0.055ng/mL
Entrez	403461
Expiration Date	Please enquire.