

## Product Datasheet

### General THB ELISA Kit (orb1201853)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | General Tetrahydrobiopterin (THB) ELISA Kit                   |
| <b>Reactivity</b>          | All   |
| <b>Range</b>               | 0.156-10ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | competitive inhibition  |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Sensitivity</b>         | 0.056ng/mL  |
| <b>Expiration Date</b>     | Please enquire.   |