

## Product Datasheet

### General PGI2 ELISA Kit (orb1201857)

<b>Description</b>	General Prostacyclin (PGI2) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	78.1-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	33.5pg/mL
<b>Expiration Date</b>	Please enquire.