

## Product Datasheet

### General PGE1 ELISA Kit (orb1201862)

<b>Description</b>	General Prostaglandin E1 (PGE1) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	62.5-4000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	15.6pg/mL
<b>Expiration Date</b>	Please enquire.