

## Product Datasheet

### General NE ELISA Kit (orb1201863)

<b>Description</b>	General Noradrenaline (NE) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	78.1-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	27.2pg/mL
<b>Expiration Date</b>	Please enquire.