

## Product Datasheet

### General GC ELISA Kit (orb1201866)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | General Glucocorticoid (GC) ELISA Kit  |
| <b>Reactivity</b>          | All  |
| <b>Range</b>               | 0.156-10ng/mL  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | competitive inhibition   |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sensitivity</b>         | 0.059ng/mL   |
| <b>Expiration Date</b>     | Please enquire.  |