

## Product Datasheet

### Human PDGFB ELISA Kit (orb1201918)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Platelet Derived Growth Factor Subunit B (PDGFB) ELISA Kit                                     |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 31.2-2000pg/mL   |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | sandwich   |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">P01127</a>   |
| <b>Sensitivity</b>         | 11.3pg/mL  |
| <b>Entrez</b>              | <a href="#">5155</a>   |
| <b>Expiration Date</b>     | Please enquire.  |