

Product Datasheet

Rat IP10 ELISA Kit (orb1201998)

Description	Rat Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48973
Sensitivity	11.5pg/mL
Entrez	245920
Expiration Date	Please enquire.