

Product Datasheet

Mouse IL19 ELISA Kit (orb1202028)

Description	Mouse Interleukin 19 (IL19) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8CJ70
Sensitivity	6.3pg/mL
Entrez	329244
Expiration Date	Please enquire.