

Product Datasheet

Mouse Bid ELISA Kit (orb1202127)

Description	Mouse BH3 Interacting Domain Death Agonist (Bid) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P70444
Sensitivity	0.112ng/mL
Entrez	12122
Expiration Date	Please enquire.