

Product Datasheet

Porcine IL8 ELISA Kit (orb1202133)

Description	Porcine Interleukin 8 (IL8) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P26894
Sensitivity	12.3pg/mL
Entrez	396880
Expiration Date	Please enquire.