

Product Datasheet

Human CSF2Rb ELISA Kit (orb1202192)

| | |
|----------------------------|--|
| Description | Human Colony Stimulating Factor 2 Receptor Beta (CSF2Rb) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates or other biological fluids. |
| Uniprot ID | P32927 |
| Sensitivity | 0.061ng/mL |
| Entrez | 1439 |
| Expiration Date | Please enquire. |