

Product Datasheet

Human GDF5 ELISA Kit (orb1202220)

Description	Human Growth Differentiation Factor 5 (GDF5) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P43026
Sensitivity	0.064ng/mL
Entrez	8200
Expiration Date	Please enquire.