

Product Datasheet

Mouse VEGFC ELISA Kit (orb1202245)

Description	Mouse Vascular Endothelial Growth Factor C (VEGFC) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P97953
Sensitivity	6.3pg/mL
Entrez	22341
Expiration Date	Please enquire.