

## Product Datasheet

### Human PDGFBB ELISA Kit (orb1202321)

<b>Description</b>	Human Platelet Derived Growth Factor BB (PDGFBB) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	13.4pg/mL
<b>Expiration Date</b>	Please enquire.