

## Product Datasheet

### Human IL1R1 ELISA Kit (orb1202336)

<b>Description</b>	Human Interleukin 1 Receptor Type I (IL1R1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P14778</b>
<b>Sensitivity</b>	0.049ng/mL
<b>Entrez</b>	<b>3554</b>
<b>Expiration Date</b>	Please enquire.