

Product Datasheet

Human BMP4 ELISA Kit (orb1202357)

Description	Human Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12644
Sensitivity	0.056ng/mL
Entrez	652
Expiration Date	Please enquire.