

Product Datasheet

Human ADIPOR2 ELISA Kit (orb1202365)

| | |
|----------------------------|--|
| Description | Human Adiponectin Receptor 2 (ADIPOR2) ELISA Kit |
| Reactivity | Human |
| Range | 0.781-50ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates or other biological fluids. |
| Uniprot ID | Q86V24 |
| Sensitivity | 0.24ng/mL |
| Entrez | 79602 |
| Expiration Date | Please enquire. |