

Product Datasheet

Human ADIPOR1 ELISA Kit (orb1202366)

Description	Human Adiponectin Receptor 1 (ADIPOR1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q96A54
Sensitivity	0.25ng/mL
Entrez	51094
Expiration Date	Please enquire.