

Product Datasheet

Rat ITaC ELISA Kit (orb1202372)

Description	Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
Reactivity	Rat
Range	46.8-3000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q7TNLO
Sensitivity	16.13pg/mL
Entrez	305236
Expiration Date	Please enquire.