

Product Datasheet

Canine TNFa ELISA Kit (orb1202421)

Description	Canine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Canine
Range	3.12-200pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P51742
Sensitivity	1.19pg/mL
Entrez	403922
Expiration Date	Please enquire.