

Product Datasheet

Human ENA78 ELISA Kit (orb1202453)

Description	Human Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42830
Sensitivity	0.063ng/mL
Entrez	6374
Expiration Date	Please enquire.