

Product Datasheet

Rat GAL1 ELISA Kit (orb1202482)

Description	Rat Galectin 1 (GAL1) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11762
Sensitivity	0.246ng/mL
Entrez	56646
Expiration Date	Please enquire.