

Product Datasheet

Human CST3 ELISA Kit (orb1202484)

Description	Human Cystatin 3 (CST3) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01034
Sensitivity	0.64ng/mL
Entrez	1471
Expiration Date	Please enquire.