

Product Datasheet

Human VEGFD ELISA Kit (orb1202561)

| | |
|----------------------------|--|
| Description | Human Vascular Endothelial Growth Factor D (VEGFD) ELISA Kit |
| Reactivity | Human |
| Range | 15.6-1000pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | 043915 |
| Sensitivity | 6.1pg/mL |
| Entrez | 2277 |
| Expiration Date | Please enquire. |