

Product Datasheet

Human IL1b ELISA Kit (orb1202591)

Description	Human Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01584
Sensitivity	12.8pg/mL
Entrez	3553
Expiration Date	Please enquire.