

Product Datasheet

Rat Gamma-Secretase Activating Protein (gSAP) ELISA Kit (orb1202659)

Description	Rat Gamma-Secretase Activating Protein (gSAP) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	D3ZPQ9
Sensitivity	0.113ng/mL
Entrez	311984
Expiration Date	Please enquire.