

## Product Datasheet

# Monkey Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit (orb1202679)

<b>Description</b>	Monkey Macrophage Inflammatory Protein 1 Beta (MIP1b) ELISA Kit
<b>Reactivity</b>	Monkey
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>Q8HYQ2</b>
<b>Sensitivity</b>	9.5pg/mL
<b>Entrez</b>	<b>574217</b>
<b>Expiration Date</b>	Please enquire.