

Product Datasheet

Monkey Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit (orb1202680)

Description	Monkey Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit
Reactivity	Monkey
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P61274
Sensitivity	0.054ng/mL
Entrez	102135739
Expiration Date	Please enquire.