

Product Datasheet

Bovine Lipopolysaccharide Binding Protein (LBP) ELISA Kit (orb1202716)

Description	Bovine Lipopolysaccharide Binding Protein (LBP) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q2TB10
Sensitivity	0.233ng/mL
Entrez	512242
Expiration Date	Please enquire.