

Product Datasheet

Bovine Cathepsin D (CTSD) ELISA Kit (orb1202728)

Description	Bovine Cathepsin D (CTSD) ELISA Kit
Reactivity	Bovine
Range	3.12-200ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1MMR6
Sensitivity	0.88ng/mL
Entrez	282883
Expiration Date	Please enquire.